

Engineering Physics Laboratory

(Course Code : PHY119)

Experiment Number: 5

AIM: Determination of the velocity of Ultrasonic waves using Ultrasonic interferometer. Hence find the compressibility of the given liquid.



Learning Objectives

- a. To understand about standing waves in liquid column.**
- b. To perform the simulated version of the actual ultrasonic interferometer experiment.**
- c. To enhance the knowledge of production of ultrasonic waves.**
- d. To understand the dependence of the velocity of ultrasonic waves on the density of the medium.**


Pre-requisite Understanding

Standing waves and resonance



- At ordinary frequencies, waves travel backwards and forwards along the string or any media/vacuum.
- Each new reflected wave has a new phase.




Standing waves and resonance



•However, at certain special frequencies, the interference produces strong standing wave patterns.




•Such a standing wave is said to be produced at resonance.





•These certain frequencies are called resonant frequencies.

Node and anti-node in a standing wave



Anti-node: Region of maximum amplitude in a standing wave

Node: Region of minimum amplitude in a standing wave



Classified of Sound Waves:

- Depends upon Frequency
- Divided into 3 groups.

Description	Frequency range Hz	Example
Infrasound	0 – 20	Earth quake
Audible	20 – 20,000	Speech, music
Ultrasonic	> 20,000 to 5M	Bat, Quartz crystal

Production of ultrasonic waves by inverse Piezo-electric effect

Principle : Inverse piezo electric effect


- If mechanical pressure is applied to one pair of opposite faces of certain crystals like quartz, equal and opposite electrical charges appear across its other faces. This is called as piezo-electric effect.
- The converse of piezo electric effect is also true.
- If an electric field is applied to one pair of faces, the corresponding changes in the dimensions of the other pair of faces of the crystal are produced. This is known as *inverse piezo electric effect* or *electrostriction*.

Step by Step guide to perform the experiment in Virtual lab

(1) Follow the slides below.

(2) Click on the YouTube video link: [http://
www.youtube.com/v/MD_zkNzF3eA&hl=en&fs=
1&rel=0&hd=1&autoplay=1](http://www.youtube.com/v/MD_zkNzF3eA&hl=en&fs=1&rel=0&hd=1&autoplay=1)

Type this link on the address bar or Click on this link : <https://www.vlab.co.in/>



The screenshot shows a web browser window with the address bar displaying www.vlab.co.in. The page content includes a header with the text "Virtual Labs" and a description: "To provide remote-access to Labs in various disciplines of Science and Engineering. These Virtual Labs would cater to students at the undergraduate level, post ...". Below this is a search bar with the text "Results from vlab.co.in" and a magnifying glass icon. The main content area is organized into two columns, each with a category header and a list of links.

www.vlab.co.in ▾

Virtual Labs

To provide remote-access to Labs in various disciplines of Science and Engineering. These **Virtual Labs** would cater to students at the undergraduate level, post ...

Results from vlab.co.in 🔍

Electronics & Communications Reference Books. DIGITAL DESIGN ? Third Edition , M ...	Physical Sciences Physical Sciences. Home » Broad Areas of Virtual Labs. Labs ...
Civil Engineering Civil Engineering. Home » Broad Areas of Virtual Labs. Labs ...	Computer Science ... Reference Books. Fundamentals of Data structures by Ellis ...
Electrical Engineering Electrical Engineering. Home » Broad Areas of Virtual Labs ...	Mechanical Engineering Mechanical Systems and Signal Processing Lab. Reference ...

This page will open. Scroll down the page end click on Physical sciences



An MHRD Govt of India Initiative

An Initiative of
Ministry of Human Resource Development
Under the National Mission on Education through **ICT**



HOME

ABOUT US

NEWSLETTER

NEW LAB DEVELOPMENT

BECOME NODAL CENTER

OUTREACH PORTAL

PARTICIPATING INSTITUTES

CONTACT US



OBJECTIVES

THE PHILOSOPHY

SALIENT FEATURES

Objectives

Scroll down the page end click on Physical sciences

Broad Areas of Virtual Labs

- ❏ Electronics & Communications
- ❏ Computer Science & Engineering
- ❏ Electrical Engineering
- ❏ Mechanical Engineering
- ❏ Chemical Engineering
- ❏ Biotechnology and Biomedical Engineering
- ❏ Civil Engineering
- ❏ Physical Sciences
- ❏ Chemical Sciences

Participating Institutes



IIT GUWAHATI



IIT DELHI



IIT BOMBAY



IIT KANPUR



IIIT
HYDERABAD



AMRITA
VISHWA
VIDYAPEETHAM



IIT MADRAS



DAYALBAGH
EDUCATIONAL
INSTITUTE

Then click on – Harmonic motion and waves lab

← → ↻ 🏠 🔒 vlab.co.in/broad-area-physical-sciences ★ ⓘ

HOME	ABOUT US	NEWSLETTER	NEW LAB DEVELOPMENT	BECOME NODAL CENTER	OUTREACH PORTAL	PARTICIPATING INSTITUTES	CONTACT US
Reference Books ▼ Syllabus Mapping ▼							
Virtual Anthropology Lab				IIT GUWAHATI			
Reference Books ▼ Syllabus Mapping ▼							
Advanced Mechanics Virtual Lab				AMRITA VISHWA VIDYAPEETHAM			
Reference Books ▼ Syllabus Mapping ▼							
Laser Optics Virtual Lab				AMRITA VISHWA VIDYAPEETHAM			
Reference Books ▼ Syllabus Mapping ▼							
Heat & Thermodynamics Lab				AMRITA VISHWA VIDYAPEETHAM			
Reference Books ▼ Syllabus Mapping ▼							
Modern Physics Lab				AMRITA VISHWA VIDYAPEETHAM			
Reference Books ▼ Syllabus Mapping ▼							
<u>Harmonic Motion and Waves Lab</u>				AMRITA VISHWA VIDYAPEETHAM			
Reference Books ▼ Syllabus Mapping ▼							

Then click on – Ultrasonic Interferometer

Not secure | vlab.amrita.edu/index.php?sub=1&brch=201



VALUE @ Amrita

Virtual Amrita Laboratories Universalizing Education



[Home](#) [Project](#) [Workshop](#) [Nodal Centres](#) [News & Events](#) [Publications](#) [Survey](#) [Contact us](#) [Login](#)

Google Custom Search

you are here->[home](#)->[physical sciences](#)->[harmonic motion and waves](#) [virtual lab](#)

Harmonic Motion and Waves Virtual Lab

Harmonic Motion and Wave lab is the interdisciplinary science that deals with the study of sound, ultrasound and infrasound (all mechanical waves in gases, liquids, and solids). The study of this lab revolves around the generation, propagation and reception of mechanical waves and vibrations.



Astable multivibrator

Design and set up an astable multivibrator and calculate the time period and duty cycle of the output wave form.



Melde's String Apparatus

Aim is to determine the frequency of electrically maintained tuning fork.



Kundt's Tube Apparatus

Kundt's tube is a device used to calculate the velocity of sound through a metal rod. Also the Young's modulus of the rod can be calculated.



Ultrasonic Interferometer

Ultrasonic Interferometer is a device used to calculate the velocity of ultrasonic sound through different liquid media.



Doppler Effect

The Doppler Effect is the perceived change in frequency of sound emitted by a source moving relative to the observer. The effect was first noted by Christian Doppler in 1842.




A.C Sonometer

Every object has a natural frequency of vibration. If kinetic energy is applied to an object at a rate that matches its natural frequency, resonance occurs and the object vibrates. In this experiment a small current, produced by a signal generator, causes

Details of the experiment will be available to you . You may login also.

- **Read the theory**
- **Procedure**
- **Complete the self evaluation to check your understanding**
- **Then click on the simulator**

Ultrasonic Interferometer



Theory



Procedure



Self Evaluation



Simulator



Assignment



Reference



Feedback

Aim:

To calculate the velocity of ultrasonic sound through different liquid media.
To calculate the adiabatic compressibility of the given liquid.

Apparatus:

Ultrasonic interferometer, sample liquids, high frequency generator etc.

Theory:


Ultrasonic interferometer is a simple device which yields accurate and consistent data, from which one can determine the velocity of ultrasonic sound in a liquid medium.

Ultrasonics:

Ultrasonic sound refers to sound pressure with a frequency greater than the human audible range (20Hz to 20 KHz). When an ultrasonic wave propagates through a medium, the molecules in that medium vibrate over very short distance in a direction parallel to the longitudinal wave. During this vibration, momentum is transferred among molecules. This causes the wave to pass through the medium.

Generation of ultrasound:

Ultrasonic can be produced by different methods. The most common methods include:




Click on the simulator

Alternative link, if adobe flash is not supported by your browser. It works in html.




http://hmw-au.vlabs.ac.in/harmonic-motion-waves/Ultrasonic_Interferometer/experiment.html

Click on the cross section icon to see the interferometer cell and switch on the frequency generator by clicking 'Power on' icon.





Ultrasonic Interferometer



Variables

Choose medium:

Water($\rho=996.458 \text{ kg/m}^3$)

Frequency of wave(MHz): 2

1 5

Reset

☒ Show cross section

☐ Show graph

☐ Show result


Result

Selection of medium and frequency

Select the medium from the 'choose medium' drop down menu to 'Water'.


Move the frequency sliding bar to 2 MHz

Selection of 'Gain' and 'ADJ' knob



Click here

Keep clicking the right arrow of the 'ADJ' knob until the ammeter reading reaches a higher value (80 in figure) from it's 'zero' position.



Click here


Keep clicking the right arrow of the 'GAIN' knob until the ammeter reading reaches a higher value (40 in figure) from it's 'zero' position. **Make sure the 'GAIN' reading is less than that of 'ADJ' reading.**

Selection of graph option



Pop-up window to observe the graph.

Click the 'Show graph' option to see the graph of the ammeter readings. A small pop up window will open.




Initially the circular scale is at the 'zero' position of the linear scale reading. You need to click the right arrow icon to move the circular scale upward.

Micro-meter scale


The micro-meter scale, as shown, consists of two parts. Linear scale and circular scale, just like any screw gauge. The circular scale moves upward by clicking the right arrow icon and it moves downward by clicking the left arrow icon, as shown the figure. The least count of this scale is 0.01 mm and all readings are in mm scale.

Tabulation of Micro-meter readings



- In the simulator, right and left arrows are provided to move the circular scale upward and downward across the linear scale respectively. Increase the micrometer setting till the anode current in the ammeter shows a new maximum. (After the first few clicks, if you click and hold the arrow, the micrometer setting will increase continuously. A single click increases it by a small increment.) Note down the micrometer reading at the new maximum.

Tabulation of Micro-meter readings



- Micro-meter reading = Linear scale reading + Circular scale reading X Least count
- Example : Linear scale reading = 1.2 mm, Circular scale reading = 35.

$$\text{Micro-meter reading} = 1.2 + (35 \times 0.01) = 1.55 \text{ mm}$$

- Note down the micrometer readings for each peak position of the current in the graph.
- These readings are the at the current readings at each anti-nodes inside the water column and distance between two anti-nodes is half of the wavelength.

Tabulation of Micro-meter readings

Order of Maxima (n)	Micrometer reading for maximum			Order of Maxima (n)	
	MSR (mm)	CSD	MSR + (CSD x LC) mm	MSR (mm)	CSD
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
Order of Maxima (n)			Micrometer reading for maximum		
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
x_{n+1} and x_n are the micro-meter readings for two consecutive maxima current values when moving the circular scale in a particular direction. And the difference between them is the half of the wavelength, multiplying that value by 2, we get the wavelength of ultrasonic wave in mm. One need to calculate the average value of the wavelength and write that down in the table.			x_{n+1} and x_n are the micro-meter readings for two consecutive maxima current values when moving the circular scale in a particular direction. And the difference between them is the half of the wavelength, multiplying that value by 2, we get the wavelength of ultrasonic wave in mm. One need to calculate the average value of the wavelength and write that down in the table.		
x_{n+1} and x_n are the micro-meter readings for two consecutive maxima current values when moving the circular scale in a particular direction. And the difference between them is the half of the wavelength, multiplying that value by 2, we get the wavelength of ultrasonic wave in mm. One need to calculate the average value of the wavelength and write that down in the table.			x_{n+1} and x_n are the micro-meter readings for two consecutive maxima current values when moving the circular scale in a particular direction. And the difference between them is the half of the wavelength, multiplying that value by 2, we get the wavelength of ultrasonic wave in mm. One need to calculate the average value of the wavelength and write that down in the table.		

Necessary Formulae

Frequency of the ultrasonic wave (f)

Wavelength of the ultrasonic wave (λ)

Velocity of the ultrasonic waves in the given liquid (v) = $\lambda \times f$

Compressibility = $1/\rho v^2$

Reading of adobe flash interface

Order of maxima (n)	MSR	C _s					
1	0.85						
2	1.2						
3	1.6						
4	2.1						
5	2.5						
6	2.8						

Velocity of sound wave: 1498 m/sec

Reading of adobe flash interface

Order of maxima (n)	MSR	CSI					
1	0.5						
2	0.8						
3	1.2						
4	1.7						
5	2.2						

Velocity of sound wave: 1680 m/sec

Reading of adobe flash interface

Order of maxima (n)	MSR	C					
1	0.4						
2	0.65						
3	1.1						
4	1.5						
5	1.9						
6	2.4						

Velocity of sound wave: 1472 m/sec

Reading of adobe flash interface

Order of maxima (n)	MSR	C					
1	0.4						
2	0.8						
3	1.2						
4	1.6						
5	2						
6	2.5						

Velocity of sound wave: 1560 m/sec

Reading of adobe html interface

Order of maxima (n)	MSR	C					
1	1.1						
2	1.4						
3	1.7						
4	2.2						
5	2.8						
6	3.1						

Velocity of sound wave: 1777 m/sec

Reading of adobe html interface

Order of maxima (n)	MSR	C					
1	1.1						
2	1.4						
3	1.7						
4	2.2						
5	2.7						
6	3.2						

Velocity of sound wave: 1808 m/sec

Reading of adobe html interface

Order of maxima (n)	MSR	C					
1	1						
2	1.4						
3	1.8						
4	2.2						
5	2.7						
6	3.2						

Velocity of sound wave: 2040 m/sec

Theory based questions

Q: 5.1 The human audible range is

- A. 20Hz - 20KHz**
- B. 15KHz- 2MHz**
- C. Above 2MHz**
- D. None of the above**

Q: 5.2 Ultrasonic waves have frequency is,

- A. Less than that of human audible range**
- B. In between human audible range**
- C. Greater than the upper human audible range**
- D. None of the above**

Q: 5.3 The velocity(v) of a wave is related to its wavelength(λ) by the relation,

A. $v = f\lambda$

B. $f = v\lambda$

C. $v = f/\lambda$

D. $\lambda = vf$

Q: 5.4 Among the following, which is the unit of adiabatic compressibility?

- A. N/m^2**
- B. m^2/N**
- C. Kg/Nm^2**
- D. None of these**

Q: 5.5 Compressibility of a medium is defined as the reciprocal of

- A. Young's modulus Rigidity modulus Density
Young's modulus**
- B. Rigidity modulus**
- C. Density**
- D. Bulk modulus**

Q: 5.6 Ultra-Sound waves are

- (A) Longitudinal**
- (B) Transverse**
- (C) Electromagnetic**
- (D) Only magnetic**

Q: 5.7 Ultra-Sound waves cannot pass through

- (A) Air**
- (B) Vacuum**
- (C) Solid**
- (D) Liquid**

Q: 5.8

The distance between an Anti-node and the next Anti-node of a transverse wave is _____.

(A) $\frac{\lambda}{4}$

$\frac{\lambda}{2}$

$\frac{\lambda}{8}$

Q: 5.9 In magnetostriction method,

- (A) Length of an object change in presence of a magnetic field.**
- (B) Color of an object changes in presence of a magnetic field.**
- (C) Resistance of a metal becomes zero in presence of a magnetic field.**
- (D) A paramagnetic material becomes a ferromagnetic material in presence of a magnetic field.**

Q: 5. 10 Which of the following materials could be used for the production of Ultrasonic wave,

- (A) Iron**
- (B) Nickel**
- (C) Quartz Crystal**
- (D) All of the above**

Q: 5.11 Speed of stationary waves is,

- (A) 1 m s^{-1}**
- (B) 2 m s^{-1}**
- (C) 3 m s^{-1}**
- (D) Zero**

Q: 5.12 In a stationary wave, nodes are at,

- (A) Fixed points**
- (B) Movable points**
- (C) There are no nodes**
- (D) Random points**

Activity based questions

A: 5.1 Which liquid you used during the Experiment?

A. Water

B. Sulphuric Acid

C. Petrol

D. Mercury

A: 5.2 If you rotate the micrometer screw anti-clock wise, which direction the circler scale will move?

A.Upward

B.Downward

C.Right side

D.Left Side

A: 5.3 If you keep the 'adj' value less than that of 'Gain' Value, can you perform the experiment?

A.Yes

B.No

A: 5.4 Which parameter does get recorded in the graph when we rotate the micrometer scale?

- A. Density of water**
- B. Frequency of the ultra-sonic wave**
- C. Ammeter reading**
- D. Voltage of the external power supply**

A: 5.5 How many different types of liquid you can choose in this experiment?

A.3

B.4

C.5

D.7

A: 5.6 What is the range of frequency of the interferometer you can change in the simulator?

A.3 MHZ-10 MHz

B.1 MHZ- 5 MHz

C.1 MHz-4 MHz

D.3 MHz-6 MHz

A: 5.7 How do you turn on the interferometer in the simulator?

A.It gets automatically turned on once run.

B.It is not needed.

C.One need to click 'Power on' icon in the simulator.

D.It has to be done by your PHY 119 teacher remotely.

A: 5.8 What is the distance between two consecutive graduation on the linear scale, as available in the simulator?

A.0.5 mm

B.1 mm

C.1.5 mm

D.2 mm

